

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-NFATC2IP



Numéro de catalogue: 11667-1-AP

Informations de base

Numéro de catalogue: 11667-1-AP	Numéro d'acquisition GenBank: BC021551	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul , Concentration: 247 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 84901	
Hôte: Lapin	Nom complet: nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 2 interacting protein	
Isotype: IgG	MW calculé 138 aa, 15 kDa	
Immunogen Catalog Number: AG2217	MW observés: 50 kDa	

Applications

Applications testées:
IHC, WB, ELISA

Spécificité de l'espèce:
Humain, souris

Informations générales

NFATC2IP, also named as NIP45, is a 419 amino acid protein, which contains one ubiquitin-like domain. TRAF1 is associated with a fraction of NFATC2IP in the cytoplasm and prevents its translocation to the nucleus. NFATC2IP in T-helper 2 (Th2) cells, regulates the magnitude of NFAT-driven transcription of a specific subset of cytokine genes, including IL3, IL4, IL5 and IL13, but not IL2.

Stockage

Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 11667-1-AP (NFATC2IP Antibody) at dilution of 1:100 incubated at 4 degree celsius over night.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 11667-1-AP (NFATC2IP Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 11667-1-AP (NFATC2IP Antibody) at dilution of 1:50 (under 40x lens).