

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-MPZL2



Numéro de catalogue: 11787-1-AP

Phare

7 Publications

Informations de base

Numéro de catalogue:
11787-1-AP

Taille:
150ul , Concentration: 450 µg/ml by
Nanodrop;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG2419

Numéro d'acquisition GenBank:
BC017774

Identification du gène (NCBI):
10205

Nom complet:
myelin protein zero-like 2

MW calculé
215 aa, 24 kDa

MW observés:
24-30 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:1000-1:4000

Applications

Applications testées:
WB, ELISA

Demandes citées:
FC, IF, IHC, WB

Spécificité de l'espèce:
Humain, rat, souris

Espèces citées:
Humain, souris

Contrôles positifs:

WB : cellules A431, cellules HaCaT

Informations générales

MPZL2 (myelin protein zero-like 2), also named as EVA (epithelial V-like antigen) or EVA1, is a type 1 transmembrane protein of the myelin PO protein family, containing an Ig-like V-type (immunoglobulin-like) domain and two potential glycosylation sites. MPZL2 is expressed in thymus epithelium and strongly downregulated by thymocyte developmental progression. Its capacity to mediate cell adhesion through a homophilic interaction and its selective regulation by T cell maturation might imply the participation of EVA in the earliest phases of thymus organogenesis.

Publications notables

Autrice	Pubmed ID	Journal	Application
QianZhi Ni	31997489	Cancer Sci	WB, IHC
Mieke Wesdorp	29961571	Am J Hum Genet	IF
Guney Bademci	29982980	Hum Genet	WB, IF

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 11787-1-AP (MPZL2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.