

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-GIMAP4



Numéro de catalogue: 11852-1-AP

Informations de base

Numéro de catalogue:

11852-1-AP

Taille:

150ul, Concentration: 600 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG2422

Numéro d'acquisition GenBank:

BC020657

Identification du gène (NCBI):

55303

Nom complet:

GTPase, IMAP family member 4

MW calculé

329 aa, 38 kDa

MW observés:

38 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

IP 0.5-4.0 µg for IP and 1:500-1:1000 for WB

Applications

Applications testées:

IP, WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules Jurkat,

IP : cellules Jurkat,

Informations générales

GIMAP4, also named as hIAN1 or HIMAP4, is a member of GIMAP family. Among GIMAP family, GIMAP4 is the only member reported to have true GTPase activity, and its transcription is found to be differentially regulated during early human CD4+ T helper (Th) lymphocyte differentiation. GIMAP4 has been previously connected mainly with T- and B-cell development and survival and T-cell apoptosis. The MW of this protein is 38 kDa, and Catalog# 11852-1-AP specially recognises the 38 kDa protein.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Jurkat cells were subjected to SDS PAGE followed by western blot with 11852-1-AP (GIMAP4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



IP Result of anti-GIMAP4 (IP:11852-1-AP, 3ug; Detection:11852-1-AP 1:800) with Jurkat cells lysate 1080ug.