

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-PHOSPHO2



Numéro de catalogue: 11869-1-AP

Informations de base

Numéro de catalogue:

11869-1-AP

Taille:

150ul, Concentration: 300 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG2451

Numéro d'acquisition GenBank:

BC022324

Identification du gène (NCBI):

493911

Nom complet:

phosphatase, orphan 2

MW calculé

241 aa, 28 kDa

MW observés:

28 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:200-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

souris

Contrôles positifs:

WB : tissu testiculaire de souris,

Informations générales

PHOSPHO2 belongs to the HAD-like hydrolase superfamily and PHOSPHO family. PHOSPHO2 is a phosphatase that has high activity toward pyridoxal 5'-phosphate (PLP). PHOSPHO2 is also active at much lower level toward pyrophosphate, phosphoethanolamine (PEA), phosphocholine (PCho), phospho-L-tyrosine, fructose-6-phosphate, p-nitrophenyl phosphate, and L-glycerophosphate. The molecular weight of PHOSPHO2 is 28 kDa.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse testis tissue were subjected to SDS PAGE followed by western blot with 11869-1-AP (PHOSPHO2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.