

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-NAV2



Numéro de catalogue: 11875-1-AP

Informations de base

Numéro de catalogue:
11875-1-AP

Taille:
150ul , Concentration: 133 µg/ml by
Bradford method using BSA as the
standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG2463

Numéro d'acquisition GenBank:
BC016054

Identification du gène (NCBI):
89797

Nom complet:
neuron navigator 2

MW calculé
268 kDa

Méthode de purification:

Purification par affinité contre
l'antigène

Applications

Applications testées:
WB,ELISA

Spécificité de l'espèce:
Humain

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HEK-293 cells were subjected to SDS PAGE followed by western blot with 11875-1-AP (NAV2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.