

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-VPREB3



Numéro de catalogue: 11915-1-AP

## Informations de base

<b>Numéro de catalogue:</b> 11915-1-AP	<b>Numéro d'acquisition GenBank:</b> BC020666	<b>Méthode de purification:</b> Purification par affinité contre l'antigène
<b>Taille:</b> 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	<b>Identification du gène (NCBI):</b> 29802	
<b>Hôte:</b> Lapin	<b>Nom complet:</b> pre-B lymphocyte 3	
<b>Isotype:</b> IgG	<b>MW calculé</b> 123 aa, 14 kDa	
<b>Immunogen Catalog Number:</b> AG2528		

## Applications

**Applications testées:**  
ELISA

**Spécificité de l'espèce:**  
Humain

## Informations générales

VPREB3 (formerly named 8HS20) is a protein that associates with mu chains in the endoplasmic reticulum of pre-B cell lines (PMID: 8491176; 22684248). VPREB3 is normally expressed by precursor B cells in bone marrow and by a subset of normal germinal center B cells in secondary lymphoid organs. It has been reported that VPREB3 expression is associated with B-cell tumors with c-MYC abnormalities (PMID: 20823132).

## Stockage

**Stockage:**  
Stocker à -20°C. Stable pendant un an après l'expédition.

**Tampon de stockage:**  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

