

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-PTMS



Numéro de catalogue: 12106-1-AP

Informations de base

Numéro de catalogue:
12106-1-AP

Taille:
150ul, Concentration: 133 µg/ml by
Bradford method using BSA as the
standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG2751

Numéro d'acquisition GenBank:
BC017025

Identification du gène (NCBI):
5763

Nom complet:
parathymosin

MW calculé
12 kDa

MW observés:
10-12 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:200-1:1000
IF 1:10-1:100

Applications

Applications testées:
IF, WB, ELISA

Spécificité de l'espèce:
Humain

Contrôles positifs:

WB : cellules U-937,

IF : cellules A549,

Informations générales

PTMS, also named as Parathymosin, is a small nuclear protein that can physically interact with glucocorticoid receptors (GR) and histone H1, respectively, and can serve as a coactivator of GR (PMID: 16150697; PMID: 15716277). It is widely distributed in mammalian tissues. Catalog# 12106-1-AP is a rabbit polyclonal antibody raised against the full-length of human PTMS. The MW of this protein is 12 kDa, and this antibody specially recognises the 12 kDa protein.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



U-937 cells were subjected to SDS PAGE followed by western blot with 12106-1-AP (PTMS antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of A549 cells using 12106-1-AP (PTMS antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).