

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-INPP4A



Numéro de catalogue: 12665-1-AP

Informations de base

Numéro de catalogue:

12665-1-AP

Taille:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG3346

Numéro d'acquisition GenBank:

BC028361

Identification du gène (NCBI):

3631

Nom complet:

inositol polyphosphate-4-phosphatase, type I, 107kDa

MW calculé

972 aa, 109 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

IP, ELISA

Spécificité de l'espèce:

Humain, rat, souris

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



IP result of anti-INPP4A (IP:12665-1-AP, 4ug;
Detection:12665-1-AP 1:300) with mouse brain
tissue lysate 3200 ug.