

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-p18 INK4c



Numéro de catalogue: 12993-1-AP

Informations de base

Numéro de catalogue:

12993-1-AP

Taille:

150ul , Concentration: 240 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG3655

Numéro d'acquisition GenBank:

BC016173

Identification du gène (NCBI):

1031

Nom complet:

cyclin-dependent kinase inhibitor 2C (p18, inhibits CDK4)

MW calculé

168 aa, 18 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

IHC, ELISA

Spécificité de l'espèce:

Humain, rat, souris

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunohistochemical analysis of paraffin-embedded human breast cancer using 12993-1-AP (p18 INK4c antibody) at dilution of 1:100 (under 10x lens).