

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-TEAD3



Numéro de catalogue: 13120-1-AP

2 Publications

## Informations de base

Numéro de catalogue:

13120-1-AP

Taille:

150ul, Concentration: 400 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG3875

Numéro d'acquisition GenBank:

BC027877

Identification du gène (NCBI):

7005

Nom complet:

TEA domain family member 3

MW calculé

435 aa, 49 kDa

MW observés:

50 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

IF 1:50-1:500

## Applications

Applications testées:

IF, WB, ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain, souris

Espèces citées:

souris

Contrôles positifs:

WB : cellules HCT 116, cellules HepG2, cellules MCF-7

IF : cellules MCF-7,

## Informations générales

TEAD3, others named ranscriptional enhancer factor 5(TEF5), a member of family of transcription factors that plays a key role in the Hippo signaling pathway, suggesting involve in organ and tumor suppression by restricting proliferation and promoting. TEAD3 contains the TEA/ATTS DNA-binding domain. Human TEAD3 has a strong mRNA expression in placenta and choriocarcinoma cells(JEG-3), but low in Hela, HepG2 and MCF-7. This is a rabbit polyclonal anti-TEAD3 antibody raised against part of C-terminal chain of human TEAD3.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Duo Zhang	25512392	J Cell Biol	WB
Yang Zhao	36993303	Res Sq	WB

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using 13120-1-AP (TEAD3 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Various lysates were subjected to SDS PAGE followed by western blot with 13120-1-AP (TEAD3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.