

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-PAX9



Numéro de catalogue: 13195-1-AP

## Informations de base

Numéro de catalogue:

13195-1-AP

Taille:

150ul

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG3897

Numéro d'acquisition GenBank:

BC001159

Identification du gène (NCBI):

5083

Nom complet:

paired box 9

MW calculé

341 aa, 36 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain, rat, souris

## Informations générales

PAX9, a member of the paired box-containing gene family, is closely related in its paired domain to PAX1. PAX9 is essential for the development of a variety of organs and skeletal elements, particularly the initial stage of odontogenesis. Mutations in PAX9 can cause selective tooth agenesis in humans. PAX9 is reported as a sensitive marker for deregulated differentiation of oesophageal keratinocytes, and it has a role in the normal differentiation process of internal stratified squamous epithelia.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

