

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ATCAY



Numéro de catalogue: 13257-2-AP

Informations de base

Numéro de catalogue:

13257-2-AP

Taille:

150ul, Concentration: 280 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG3658

Numéro d'acquisition GenBank:

BC026217

Identification du gène (NCBI):

85300

Nom complet:

ataxia, cerebellar, Cayman type

MW calculé

371 aa, 42 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13257-2-AP (ATCAY antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.