

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-SH3RF2



Numéro de catalogue: 13306-1-AP

Informations de base

Numéro de catalogue:	Numéro d'acquisition GenBank:	Méthode de purification:
13306-1-AP	BC031650	Purification par affinité contre l'antigène
Taille:	Identification du gène (NCBI):	Dilutions recommandées:
150ul, Concentration: 600 µg/ml by Nanodrop;	153769	WB 1:500-1:1000
Hôte:	Nom complet:	IHC 1:50-1:500
Lapin	SH3 domain containing ring finger 2	IF 1:50-1:500
Isotype:	MW calculé	
IgG	729 aa, 79 kDa	
Immunogen Catalog Number:		
AG4127		

Applications

Applications testées:	Contrôles positifs:
IF, IHC, WB, ELISA	WB : tissu cardiaque de souris, tissu cardiaque de rat
Spécificité de l'espèce:	IHC : tissu de cancer du côlon humain,
Humain, rat, souris	IF : cellules Caco-2,
Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9,0; (*) À défaut, 'le démasquage de l'antigène peut être 'effectué avec un tampon citrate pH 6,0.	

Informations générales

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




Various lysates were subjected to SDS PAGE followed by western blot with 13306-1-AP (SH3RF2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 13306-1-AP (SH3RF2 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 13306-1-AP (SH3RF2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed Caco-2 cells using SH3RF2 antibody (13306-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).