

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-TULA; UBASH3A



Numéro de catalogue: 13357-1-AP

Informations de base

Numéro de catalogue:

13357-1-AP

Taille:

150ul, Concentration: 180 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG4170

Numéro d'acquisition GenBank:

BC028138

Identification du gène (NCBI):

53347

Nom complet:

ubiquitin associated and SH3 domain containing, A

MW calculé

74 kDa, 70 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain, rat, souris

Informations générales

UBASH3A gene encodes ubiquitin-associated and SH3 domain-containing protein A (UBASH3A)/ Suppressor of T-cell receptor signaling 2 (STS2), also known as T-cell ubiquitin ligand 1 (TULA1) which is one of two family members belonging to the T-cell ubiquitin ligand (TULA) family. UBASH3A gene expression is limited to only a few tissues, with its highest expression in spleen, peripheral blood leukocytes, and bone marrow. UBASH3A protein is associated with c-Cbl and ubiquitylated proteins and has been implicated in the regulation of signaling mediated by protein-tyrosine kinases. Other family members can negatively regulate T-cell signaling. UBASH3A can facilitate growth factor withdrawal-induced apoptosis in T cells via its interaction with apoptosis-inducing factor (AIF).

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

