

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-ANGPTL5



Numéro de catalogue: 13986-1-AP

## Informations de base

Numéro de catalogue:  
13986-1-AP

Taille:  
150ul, Concentration: 500 µg/ml by  
Nanodrop and 340 µg/ml by Bradford  
method using BSA as the standard;

Hôte:  
Lapin

Isotype:  
IgG

Immunogen Catalog Number:  
AG5053

Numéro d'acquisition GenBank:  
BC049170

Identification du gène (NCBI):  
253935

Nom complet:  
angiopoietin-like 5

MW calculé  
44 kDa

MW observés:  
45-47 kDa

Méthode de purification:  
Purification par affinité contre  
l'antigène

Dilutions recommandées:  
WB 1:500-1:2000

## Applications

Applications testées:  
WB, ELISA

Spécificité de l'espèce:  
Humain

Contrôles positifs:

WB : cellules HEK-293,

## Informations générales

Angiopoietin-like protein 5 (ANGPTL5) is a member of the angiopoietin-like proteins (ANGPTLs) family. Eight members have been described so far (ANGPTL1 through ANGPTL8), and have been shown to have various physiological functions in lipid and glucose metabolism, inflammation, hematopoiesis, angiogenesis and cancer.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



HEK-293 cells were subjected to SDS PAGE followed by western blot with 13986-1-AP (ANGPTL5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.