

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-NR1H2



Numéro de catalogue: 14278-1-AP

Phare

9 Publications

## Informations de base

Numéro de catalogue: 14278-1-AP	Numéro d'acquisition GenBank: BC047750	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul, Concentration: 550 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 7376	
Hôte: Lapin	Nom complet: nuclear receptor subfamily 1, group H, member 2	
Isotype: IgG	MW calculé 51 kDa	
Immunogen Catalog Number: AG5635	MW observés: 50 kDa, 40 kDa	

## Applications

Applications testées:  
ELISA

Demandes citées:  
IF, IP, WB

Spécificité de l'espèce:  
Humain, rat, souris

Espèces citées:  
Humain, souris

## Informations générales

NR1H2, also known as LXRB, is a liver X receptors, which form a subfamily of the nuclear receptor superfamily and are key regulators of macrophage function, controlling transcriptional programs involved in lipid homeostasis and inflammation. It preferentially bind to double-stranded oligonucleotide direct repeats, which have the consensus half-site sequence 5'-AGGTCA-3' and 4-nt spacing (DR-4). It involves in the regulation of cholesterol uptake through MYLIP-dependent ubiquitination of LDLR, VLDLR and LRP8.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Tain-Junn Cheng	27622732	Chem Biol Interact	WB
Xingzhe Ma	24947959	Int J Cancer	WB
Qixue Wang	27250582	Sci Rep	WB

## Stockage

Stockage:  
Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

