

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-ASK1



Numéro de catalogue: 14385-1-AP

3 Publications

## Informations de base

Numéro de catalogue:  
14385-1-AP

Taille:  
150ul, Concentration: 133 µg/ml by  
Bradford method using BSA as the  
standard;

Hôte:  
Lapin

Isotype:  
IgG

Immunogen Catalog Number:  
AG5818

Numéro d'acquisition GenBank:  
BC054503

Identification du gène (NCBI):  
4217

Nom complet:  
mitogen-activated protein kinase  
kinase kinase 5

MW calculé  
155 kDa

Méthode de purification:  
Purification par affinité contre  
l'antigène

## Applications

Applications testées:  
ELISA

Demandes citées:  
WB

Spécificité de l'espèce:  
Humain, rat, souris

Espèces citées:  
Humain, rat, singe

## Informations générales

ASK1 (Apoptosis signal-regulating kinase 1) is also named as MAP3K5, MAPKKK5, MEKK5 and belongs to the MAP kinase kinase kinase subfamily. It is an evolutionarily conserved mitogen activated protein 3-kinase that activates both Jnk and p38 mitogen-activated protein kinases. ASK1 is activated in response to various cytotoxic stresses including TNF, Fas and reactive oxygen species (ROS) such as H<sub>2</sub>O<sub>2</sub>, and activates c-Jun NH 2-terminal kinase (JNK) and p38.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Xiaoye Liu	27762484	Cell Microbiol	WB
Na Zhao	29522763	J Mol Cell Cardiol	WB
Chen Xiaoping X	20619274	Life Sci	WB

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

