

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ITM2A



Numéro de catalogue: 14407-1-AP

Phare

2 Publications

Informations de base

Numéro de catalogue:

14407-1-AP

Taille:

150ul, Concentration: 400 µg/ml by Nanodrop and 353 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG5692

Numéro d'acquisition GenBank:

BC040437

Identification du gène (NCBI):

9452

Nom complet:

integral membrane protein 2A

MW calculé

30 kDa

MW observés:

30 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

Applications

Applications testées:

WB, ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

souris

Contrôles positifs:

WB : tissu pulmonaire de souris,

Informations générales

Integral membrane protein 2A (ITM2A) is a single-pass type II membrane protein belonging to a family, which is composed of at least two other members ITM2B and ITM2C. ITM2A can be induced during thymocyte selection and T cell activation. In early stages of differentiation ITM2A may be associated with an inhibition of the initiation of chondrogenesis. Catalog#14407-1-AP recognises two band of 30 kDa and 45-48 kDa, and the additional 45-48 kDa band may be due to glycosylation.

Publications notables

Autrice	Pubmed ID	Journal	Application
Rui Zhang	33680917	Front Oncol	WB
Stephanie J Davies	28872940	Adipocyte	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse lung tissue were subjected to SDS PAGE followed by western blot with 14407-1-AP (ITM2A antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.