

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-SNX11



Numéro de catalogue: 14953-1-AP

Informations de base

Numéro de catalogue:
14953-1-AP

Taille:
150ul, Concentration: 600 µg/ml by
Nanodrop and 287 µg/ml by Bradford
method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG6791

Numéro d'acquisition GenBank:
BC000768

Identification du gène (NCBI):
29916

Nom complet:
sorting nexin 11

MW calculé
30 kDa

MW observés:
35 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:200-1:1000

Applications

Applications testées:
WB, ELISA

Spécificité de l'espèce:
Humain, rat, souris

Contrôles positifs:

WB: cellules HeLa,

Informations générales

SNX11 belongs to the sorting nexin family of proteins which contain a phox (PX) domain. PX domain is a phosphoinositide binding domain. SNX family members are involved in intracellular trafficking. SNX11 protein lacks C-terminal coiled-coil motifs found in several other SNX family members.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HeLa cells were subjected to SDS PAGE followed by western blot with 14953-1-AP (SNX11 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.