

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-HDAC7



Numéro de catalogue: 15322-1-AP

1 Publications

Informations de base

Numéro de catalogue:

15322-1-AP

Taille:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG1101

Numéro d'acquisition GenBank:

BC006453

Identification du gène (NCBI):

51564

Nom complet:

histone deacetylase 7

MW calculé

103 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

rat

Informations générales

Histone deacetylases (HDAC) are a class of enzymes that remove the acetyl groups from the lysine residues leading to the formation of a condensed and transcriptionally silenced chromatin. At least 4 classes of HDAC were identified. Just like other class II HDACs, HDAC7 shuttles in and out of the nucleus depending on different signals. This gene encoding human HDAC7 is orthologous to mouse HDAC7 gene whose protein promotes repression mediated via the transcriptional corepressor SMRT. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. This antibody is a rabbit polyclonal antibody raised against residues near the C terminus of human HDAC7.

Publications notables

Autrice	Pubmed ID	Journal	Application
He Min M	23480850	J Cell Mol Med	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

