

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-HEYL



Numéro de catalogue: 15679-1-AP

1 Publications

Informations de base

Numéro de catalogue:

15679-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG8163

Numéro d'acquisition GenBank:

BC006087

Identification du gène (NCBI):

26508

Nom complet:

hairy/enhancer-of-split related with YRPW motif-like

MW calculé

328 aa, 35 kDa

MW observés:

35-40 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

IP 0.5-4.0 ug for IP and 1:200-1:1000 for WB

IF 1:10-1:100

Applications

Applications testées:

IF, IP, WB, ELISA

Demandes citées:

IF

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain

Contrôles positifs:

WB : cellules MDA-MB-453s, cellules HeLa

IP : cellules HeLa,

IF : cellules HeLa,

Informations générales

HEYL, also named as Hairy/enhancer-of-split related with YRPW motif-like protein, is a 328 amino acid protein, which belongs to the HEY family. HEYL as downstream effector of Notch signaling may be required for cardiovascular development. HEYL represses transcription by the cardiac transcriptional activators GATA4 and GATA6.

Publications notables

Autrice	Pubmed ID	Journal	Application
Hao Sun	37026620	Adv Sci (Weinh)	IF

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées



Immunofluorescent analysis of HeLa cells, using HEYL antibody 15679-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with 15679-1-AP (HEYL Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



IP Result of anti-HEYL (IP:15679-1-AP, 4ug; Detection:15679-1-AP 1:300) with HeLa cells lysate 2600ug.