

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-GRB7



Numéro de catalogue: 15738-1-AP

## Informations de base

Numéro de catalogue:

15738-1-AP

Taille:

150ul , Concentration: 167 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG8537

Numéro d'acquisition GenBank:

BC006535

Identification du gène (NCBI):

2886

Nom complet:

growth factor receptor-bound protein

7

MW calculé

532 aa, 60 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain, rat, souris

## Informations générales

Grb7 belongs to a family of adapter proteins. This family, including Grb10, Grb7 and Grb 14, contains a proline-rich, a RA, a PH, a BPS and SH2 domain. Grb7 associate with a number of receptors including EGFR and FAK. Functions of Grb7 include mobility because it associates receptors such FAK, and direct evidence from overexpression of this proteins. This antibody recognize the C-terminal of GRB7.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

