

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-LOC81691



Numéro de catalogue: 15762-1-AP

Informations de base

Numéro de catalogue:

15762-1-AP

Taille:

150ul, Concentration: 193 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG8393

Numéro d'acquisition GenBank:

BC007646

Identification du gène (NCBI):

81691

Nom complet:

exonuclease NEF-sp

MW calculé

774 aa, 87 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



human kidney tissue were subjected to SDS PAGE followed by western blot with 15762-1-AP (LOC81691 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.