

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-FRK



Numéro de catalogue: 15936-1-AP

## Informations de base

Numéro de catalogue:

15936-1-AP

Taille:

150ul , Concentration: 340 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG10024

Numéro d'acquisition GenBank:

BC012916

Identification du gène (NCBI):

2444

Nom complet:

fyn-related kinase

MW calculé

505 aa, 56 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain, rat, souris

## Informations générales

The fyn-related kinase (FRK) belongs to the tyrosine kinase family of protein kinases and is also named as PTK5 (protein-tyrosine kinase 5), RAK (Nuclear tyrosine protein kinase). The Frk is a tyrosine kinase expressed in various tissues, such as the gastrointestinal tract and pancreatic islet cells. It affects pancreatic beta cell number during embryogenesis and promotes beta cell cytotoxic signals in response to streptozotocin (PMID:20157399). It may function as a tumor suppressor (PMID:19345329).

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

