

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-FIP200



Numéro de catalogue: 16172-1-AP

1 Publications

## Informations de base

Numéro de catalogue:

16172-1-AP

Numéro d'acquisition GenBank:

NM\_001083617

Méthode de purification:

Purification par affinité contre l'antigène

Taille:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Identification du gène (NCBI):

9821

Nom complet:

RB1-inducible coiled-coil 1

Hôte:

Lapin

MW calculé

183 kDa

Isotype:

IgG

MW observés:

200 kDa

## Applications

Applications testées:

ELISA

Demandes citées:

IHC

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain

## Informations générales

RB1CC1, also named as RBICC, is implicated in the regulation of RB1 expression and functions as a DNA-binding transcription factor. It is a potent regulator of the RB1 pathway and a mediator that plays a crucial role in muscular differentiation. Its expression is, thus, a prerequisite for myogenic differentiation. Involved in autophagy. RB1CC1 is required for autophagosome formation. It is probably involved in the tumorigenesis of breast cancer. RB1CC1 is frequently mutated in breast cancer and shows characteristics of a classical tumor suppressor gene. This antibody is a rabbit polyclonal antibody raised against a 18Aa peptide mapping within human RB1CC1.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Lei Wang	25689718	Biochem Biophys Res Commun	IHC

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

