

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-STK32A



Numéro de catalogue: 16223-1-AP

Informations de base

Numéro de catalogue:

16223-1-AP

Taille:

150ul, Concentration: 600 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG9276

Numéro d'acquisition GenBank:

BC021666

Identification du gène (NCBI):

202374

Nom complet:

serine/threonine kinase 32A

MW calculé

396aa, 46 kDa; 166aa, 20 kDa

MW observés:

46 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, rat, souris

Contrôles positifs:

WB: cellules HeLa, cellules MCF-7

Informations générales

STK32A (Serine/threonine-protein kinase 32A) is also named as YANK1 and belongs to Ser/Thr protein kinase family. It can affect cell metabolism by protein serine/threonine kinase phosphorylation. The gene encoding STK32A maps to human chromosome 5 and this gene is associated with Cockayne syndrome through the ERCC8 gene and familial adenomatous polyposis through the adenomatous polyposis coli (APC) tumor suppressor gene.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HeLa cells were subjected to SDS PAGE followed by western blot with 16223-1-AP (STK32A antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.