

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-SLC17A3



Numéro de catalogue: 16404-1-AP

## Informations de base

Numéro de catalogue:

16404-1-AP

Taille:

150ul, Concentration: 93 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG9526

Numéro d'acquisition GenBank:

BC017952

Identification du gène (NCBI):

10786

Nom complet:

solute carrier family 17 (sodium phosphate), member 3

MW calculé

420 aa, 46 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain

## Informations générales

SLC17A3 is one member of type I Na/Phosphate co-transporter, also named as NPT4. Rat SLC17A3 was reported as 55/46 kD protein in liver and kidney tissues. The two bands were from different transcripts of the gene. This antibody may recognize a 66kd band, which may be SLC17A1(NPT1), a homologous of SLC17A3.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

