

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZNHIT2



Numéro de catalogue: 16885-1-AP

Informations de base

Numéro de catalogue:
16885-1-AP

Taille:
150ul, Concentration: 450 µg/ml by
Nanodrop and 220 µg/ml by Bradford
method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG10351

Numéro d'acquisition GenBank:
BC052240

Identification du gène (NCBI):
741

Nom complet:
zinc finger, HIT type 2

MW calculé
403 aa, 43 kDa

MW observés:
45-50 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:500-1:2000

Applications

Applications testées:
WB, ELISA

Spécificité de l'espèce:
Humain, souris

Contrôles positifs:

WB : cellules HeLa, cellules K-562, tissu cérébral de
souris, tissu cérébral humain

Informations générales

ZNHIT2 belongs to a family of HIT-finger proteins characterized by the presence of a HIT domain, a roughly 50-residue zinc-binding domain containing conserved cysteine (Cys) and histidine (His) residues. ZNHIT2 is a cofactor of R2TP/PFDL that targets the U5 snRNP. ZNHIT2 is highly expressed in testis.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 16885-1-AP (ZNHIT2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.