

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZNF22



Numéro de catalogue: 17063-1-AP

1 Publications

Informations de base

Numéro de catalogue:
17063-1-AP

Taille:
150ul, Concentration: 350 µg/ml by Nanodrop and 260 µg/ml by Bradford method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG10833

Numéro d'acquisition GenBank:
BC010642

Identification du gène (NCBI):
7570

Nom complet:
zinc finger protein 22 (KOX 15)

MW calculé:
224 aa, 26 kDa

MW observés:
26 kDa

Méthode de purification:
Purification par affinité contre l'antigène

Dilutions recommandées:
WB 1:200-1:1000

Applications

Applications testées:
WB, ELISA

Demandes citées:
WB

Spécificité de l'espèce:
Humain, rat, souris

Espèces citées:
Humain

Contrôles positifs:

WB: cellules Jurkat,

Informations générales

ZNF22, also named as Zinc finger protein 22, is a 224 amino acid protein, which contains 5 C2H2-type zinc fingers and belongs to the krueppel C2H2-type zinc-finger protein family. In the embryo, ZNF22 is expressed in developing craniofacial structures including dental epithelium of maxillary molar tooth organs, tongue epithelium and muscle, and craniofacial bone osteoblasts. In the adult, ZNF22 is expressed in mesoderm-derived tissues such as skeletal muscle, heart, kidney and liver. Intermediate expression in spleen, thymus and brain. ZNF22 binds DNA through the consensus sequence 5'-CAATG-3' and may be involved in transcriptional regulation and may play a role in tooth formation.

Publications notables

Autrice	Pubmed ID	Journal	Application
Baicheng Zhu	36518188	Front Mol Neurosci	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Jurkat cells were subjected to SDS PAGE followed by western blot with 17063-1-AP (ZNF22 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.