

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-TCP11L1



Numéro de catalogue: 17493-1-AP

Informations de base

Numéro de catalogue: 17493-1-AP	Numéro d'acquisition GenBank: BC041696	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul, Concentration: 350 µg/ml by Nanodrop and 253 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 55346	Dilutions recommandées: WB 1:500-1:1000
Hôte: Lapin	Nom complet: t-complex 11 (mouse)-like 1	
Isotype: IgG	MW calculé: 509 aa, 57 kDa	
Immunogen Catalog Number: AG11612	MW observés: 57 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules Jurkat,
Spécificité de l'espèce: Humain	

Informations générales

T-complex protein 11-like protein 1 (TCP11L1) is a member of the TCP11 family. It has been reported to share its functional domain with and has similar subcellular localization to the TCP11 protein (PMID: 33568072; Liu Y et al., Chin J Appl Environ Biol. 2008;14:783-6.). A study has described a patient with cryptorchidism and nonobstructive azoospermia presenting a novel microdeletion at 11p13 in TCP11L1 gene (PMID: 24912414).

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Jurkat cells were subjected to SDS PAGE followed by western blot with 17493-1-AP (TCP11L1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.