

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-NOL11



Numéro de catalogue: 17613-1-AP

Informations de base

Numéro de catalogue: 17613-1-AP	Numéro d'acquisition GenBank: BC064404	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul , Concentration: 400 µg/ml by Nanodrop;	Identification du gène (NCBI): 25926	Dilutions recommandées: WB 1:500-1:1000
Hôte: Lapin	Nom complet: nucleolar protein 11	
Isotype: IgG	MW calculé: 719 aa, 81 kDa	
Immunogen Catalog Number: AG11783	MW observés: 75-81 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : tissu cérébral de souris, tissu cérébral de rat
Spécificité de l'espèce: Humain, rat, souris	

Informations générales

NOL11, also named as L14, is a ribosome biogenesis factor. It is required for both optimal rDNA transcription and small subunit (SSU) pre-rRNA processing at sites A', A0, 1 and 2b. NOL11 forms a protein complex with two tryptophan-aspartic acid (WD) repeat proteins named WD-repeat protein 43 (WDR43) and Cirhin in mitotic cells. NOL11 plays a role in pathogenesis of NAIC. There're some isoforms with MW 81 kDa and 19 kDa.

Stockage

Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 17613-1-AP (NOL11 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.