

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-Thrombospondin 1



Numéro de catalogue: 18304-1-AP

Phare

28 Publications

Informations de base

Numéro de catalogue:

18304-1-AP

Taille:

150ul, Concentration: 600 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

BC028145

Identification du gène (NCBI):

7057

Nom complet:

thrombospondin 1

MW calculé

129 kDa

MW observés:

180 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

IHC 1:50-1:500

Applications

Applications testées:

IHC, WB, ELISA

Demandes citées:

IF, IHC, WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain, rat, souris

Contrôles positifs:

WB : cellules C6, cellules A431, cellules HeLa, cellules NIH/3T3, plaquettes de sang périphérique humain

IHC : tissu de cancer du poumon humain, tissu hépatique humain, tissu placentaire humain

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9,0; (*) À défaut, le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6,0.

Informations générales

Thrombospondin 1 (TSP-1), an extracellular matrix protein, is the first identified natural angiogenesis inhibitor. A variety of normal cells, including endothelial cells, fibroblasts, adipocytes, smooth muscle cells, monocytes, macrophages, and transformed cells such as malignant glioma cells, secrete Thrombospondin 1. It regulates cellular phenotype during tissue genesis and repair. It acts as a molecular facilitator by bringing together cytokines, growth factors, matrix components, membrane receptors and extracellular proteases. Thrombospondin 1 binds to a wide variety of integrin and non-integrin cell surface receptors. It is also a major activator of transforming growth factor (TGFβ1). Secreted TSP1 is a glycoprotein with a molecular mass of 150-180 kDa.

Publications notables

Autrice	Pubmed ID	Journal	Application
Gang Xu	34745323	Comput Math Methods Med	WB
Liming Zhang	34786057	Am J Transl Res	WB,IHC
Kun Jia	26993605	Oncotarget	IHC

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 18304-1-AP (Thrombospondin 1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 18304-1-AP (Thrombospondin 1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 18304-1-AP (Thrombospondin 1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).