

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-HOXD11



Numéro de catalogue: 18734-1-AP

3 Publications

## Informations de base

Numéro de catalogue:

18734-1-AP

Taille:

150ul, Concentration: 450 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

BC160135

Identification du gène (NCBI):

3237

Nom complet:

homeobox D11

MW calculé

35 kDa

MW observés:

32-36 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:200-1:1000

## Applications

Applications testées:

WB, ELISA

Demandes citées:

IF, IHC, WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain

Contrôles positifs:

WB : cellules SMMC-7721,

## Informations générales

HOXD11, also named as HOX4F, belongs to the Abd-B homeobox family. HOXD11 is a sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. HOXD11 plays a role in forelimb morphogenesis.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Karim Fekir	30837223	Cancer Res	IHC
Jian Hui Low	31303547	Cell Stem Cell	IF
Mandi Mu	37314930	Cell Rep	WB

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 18734-1-AP (HOXD11 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.