

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-SMARCA2-Specific



Numéro de catalogue: 19718-1-AP

## Informations de base

Numéro de catalogue:

19718-1-AP

Taille:

150ul

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_003070

Identification du gène (NCBI):

6595

Nom complet:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2

MW calculé

181 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain

## Informations générales

SMARCA2 is a transcriptional coactivator cooperating with nuclear hormone receptors to potentiate transcriptional activation. SMARCA2 is involved in vitamin D-coupled transcription regulation via its association with the WINAC complex, a chromatin-remodeling complex recruited by vitamin D receptor (VDR), which is required for the ligand-bound VDR-mediated transrepression of the CYP27B1 gene. The antibody is specific to SMARCA2

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

