

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-SNTB1-Specific



Numéro de catalogue: 19719-1-AP

## Informations de base

Numéro de catalogue:

19719-1-AP

Taille:

150ul, Concentration: 287 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_021021

Identification du gène (NCBI):

6641

Nom complet:

syntrophin, beta 1 (dystrophin-associated protein A1, 59kDa, basic component 1)

MW calculé

58 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain, rat, souris

## Informations générales

SNTB1, also named as SNT2B1, BSYN2 and TIP-43, belongs to the syntrophin family. SNTB1 is an adapter protein that binds to and probably organizes the subcellular localization of a variety of membrane proteins. It may link various receptors to the actin cytoskeleton and the dystrophin glycoprotein complex. This antibody is specific to SNTB1.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

