

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-SH2B2



Numéro de catalogue: 19746-1-AP

## Informations de base

Numéro de catalogue:

19746-1-AP

Taille:

150ul

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_020979

Identification du gène (NCBI):

10603

Nom complet:

SH2B adaptor protein 2

MW calculé

68 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain

## Informations générales

SH2B2, also named as APS, belongs to the SH2B adapter family. SH2B2 is an adapter protein for several members of the tyrosine kinase receptor family. It is involved in multiple signaling pathways. SH2B2 may be involved in coupling from immunoreceptor to Ras signaling. SH2B2 acts as a negative regulator of cytokine signaling in collaboration with CBL. SH2B2 binds to EPOR and suppresses EPO-induced STAT5 activation, possibly through a masking effect on STAT5 docking sites in EPOR. It suppresses PDGF-induced mitogenesis. SH2B2 may induce cytoskeletal reorganization via interaction with VAV3.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

