

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-TMEM176B



Numéro de catalogue: 19825-1-AP

2 Publications

Informations de base

Numéro de catalogue:
19825-1-AP

Taille:
150ul, Concentration: 400 µg/ml by Nanodrop and 200 µg/ml by Bradford method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG13918

Numéro d'acquisition GenBank:
BC004358

Identification du gène (NCBI):
28959

Nom complet:
transmembrane protein 176B

MW calculé
270 aa, 29 kDa

MW observés:
31 kDa, 23 kDa

Méthode de purification:
Purification par affinité contre l'antigène

Dilutions recommandées:
WB 1:200-1:1000

Applications

Applications testées:
WB, ELISA

Demandes citées:
WB

Spécificité de l'espèce:
Humain

Espèces citées:
Humain

Contrôles positifs:

WB : cellules A549, tissu placentaire humain

Informations générales

TMEM176B (also known as LR8) is a transmembrane protein belongs to the MS4A family of proteins. It may play a role in the process of maturation of dendritic cells. TMEM176B was first discovered in human lung fibroblasts and is associated with human small cell lung carcinoma. Studies suggest that TMEM176B might be involved in cancer and might serve as targets for antiangiogenic therapy of cancer patients (PMID: 24602018, 22244448). There are two isoforms produced by alternative splicing.

Publications notables

Autrice	Pubmed ID	Journal	Application
Yingying Zeng	35734424	Front Genet	WB
Mercedes Segovia	31085177	Cancer Cell	

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



A549 cells were subjected to SDS PAGE followed by western blot with 19825-1-AP (TMEM176B Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.