

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-CD133-2,3,4,5



Numéro de catalogue: 19947-1-AP

## Informations de base

Numéro de catalogue:

19947-1-AP

Taille:

150ul

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_001145847

Identification du gène (NCBI):

8842

Nom complet:

prominin 1

MW calculé

97 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain

## Informations générales

CD133, also known as PROM1 (prominin-1) or AC133, belongs to the prominin family. CD133 is a transmembrane glycoprotein with an NH<sub>2</sub>-terminal extracellular domain, five transmembrane loops and a cytoplasmic tail. The expression of CD133 has been reported in hematopoietic stem cells, endothelial progenitor cells, neuronal and glial stem cells, suggesting the potential role of CD133 as a cell surface marker of adult stem cells. CD133 has also been reported as a marker of cancer stem cells in various human tumors. This antibody recognizes transcript variant 2,3,4,5 of CD133.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

