

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-HTR2C



Numéro de catalogue: 19959-1-AP

## Informations de base

Numéro de catalogue:

19959-1-AP

Taille:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_000868

Identification du gène (NCBI):

3358

Nom complet:

5-hydroxytryptamine (serotonin) receptor 2C

MW calculé

52 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain, souris

## Informations générales

HTR2C, also named as HTR1C, belongs to the G-protein coupled receptor 1 family. HTR2C is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. It mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system. It is glycosylated at N-term. The antibody recognizes the C-term of HTR2C.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

