

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-CPXM2-Specific



Numéro de catalogue: 20082-1-AP

Informations de base

Numéro de catalogue:

20082-1-AP

Taille:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM_198148

Identification du gène (NCBI):

119587

Nom complet:

carboxypeptidase X (M14 family), member 2

MW calculé

86 kDa

MW observés:

86 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, souris

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 20082-1-AP (CPXM2-Specific Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.