

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-FOXH1-Specific



Numéro de catalogue: 20268-1-AP

Informations de base

Numéro de catalogue:
20268-1-AP

Taille:
150ul, Concentration: 160 µg/ml by
Nanodrop and 133 µg/ml by Bradford
method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Numéro d'acquisition GenBank:
NM_003923

Identification du gène (NCBI):
8928

Nom complet:
forkhead box H1

MW calculé
39 kDa

MW observés:
40-44 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:200-1:1000

Applications

Applications testées:
WB, ELISA

Spécificité de l'espèce:
Humain, rat, souris

Contrôles positifs:

WB : tissu hépatique de souris, cellules L02

Informations générales

FOXH1, also named as FAST1 and FAST2, is a transcriptional activator. FOXH1 recognizes and binds to the DNA sequence 5'-TGT[GT][GT]ATT-3'. It is required for induction of the gooseoid (GSC) promoter by TGF-beta or activin signaling. It forms a transcriptionally active complex containing FOXH1/SMAD2/SMAD4 on a site on the GSC promoter called TARE (TGF-beta/activin response element). The antibody is specific to FOXH1.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse liver tissue were subjected to SDS PAGE followed by western blot with 20268-1-AP (FOXP1-Specific antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.