

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-Factor XIIIb



Numéro de catalogue: 20269-1-AP

Informations de base

Numéro de catalogue:

20269-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM_001994

Identification du gène (NCBI):

2165

Nom complet:

coagulation factor XIII, B polypeptide

MW calculé

76 kDa

MW observés:

76-80 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:3000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : plasma humain,

Informations générales

Coagulation factor XIII (FXIII) is a blood proenzyme that catalyzes intermolecular cross-linking of fibrinogen and is involved in the clotting cascade. Factor XIII is composed of two catalytic FXIII-A and two protective/regulatory FXIII-B subunits. FXIII-B is the carrier/regulatory protein of its partner FXIII-A subunit.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



human plasma were subjected to SDS PAGE followed by western blot with 20269-1-AP (Factor XIIIb antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.