

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-TNFAIP8L3



Numéro de catalogue: 20693-1-AP

1 Publications

## Informations de base

Numéro de catalogue:

20693-1-AP

Taille:

150ul, Concentration: 600 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM\_207381

Identification du gène (NCBI):

388121

Nom complet:

tumor necrosis factor, alpha-induced protein 8-like 3

MW calculé

33 kDa

MW observés:

37 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

## Applications

Applications testées:

WB, ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain, souris

Espèces citées:

Humain

Contrôles positifs:

WB : cellules A549, cellules Jurkat, cellules MCF-7, tissu splénique de souris

## Informations générales

TNFAIP8L3 (tumor necrosis factor alpha-induced protein 8-like protein 3), also known as TIPE3, is a member of the TNFAIP8 family which also include TNFAIP8 (TIPE), TNFAIP8L1 (TIPE1) and TNFAIP8L2 (TIPE2) (PMID: 18455983; 21924498). TNFAIP8L3 is a transfer protein for phosphoinositide second messengers. It is markedly upregulated in human cancer cells and may represent a therapeutic target for treating malignant diseases (PMID: 25242044). This antibody is specific to TNFAIP8L3.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Xiaoping Huang	30084972	Hum Mol Genet	WB

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 20693-1-AP (TNFAIP8L3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.