

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-RNF212



Numéro de catalogue: 21178-1-AP

Informations de base

Numéro de catalogue:

21178-1-AP

Taille:

150ul , Concentration: 300 µg/ml by Nanodrop and 200 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG15457

Numéro d'acquisition GenBank:

BC050356

Identification du gène (NCBI):

285498

Nom complet:

ring finger protein 212

MW calculé

297 aa, 33 kDa

MW observés:

22 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, rat, souris

Contrôles positifs:

WB : tissu cardiaque humain,

Informations générales

RNF212(RING finger protein 212)acts as a regulator of crossing-over during meiosis: required to couple chromosome synapsis to the formation of crossover-specific recombination complexes.It's may be a SUMO E3 ligase that mediate sumoylation of target proteins MSH4 and MSH5, leading to enhance their binding to recombination sites.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



human heart tissue were subjected to SDS PAGE followed by western blot with 21178-1-AP (RNF212 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.