

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-PAX4



Numéro de catalogue: 21387-1-AP

Informations de base

Numéro de catalogue:

21387-1-AP

Numéro d'acquisition GenBank:

BC074761

Méthode de purification:

Purification par affinité contre l'antigène

Taille:

150ul, Concentration: 187 µg/ml by Bradford method using BSA as the standard;

Identification du gène (NCBI):

5078

Nom complet:

paired box 4

Hôte:

Lapin

MW calculé

343 aa, 37 kDa

Isotype:

IgG

Immunogen Catalog Number:

AG16022

Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain

Informations générales

PAX4, a paired-box transcription factor, is a key regulator of pancreatic islet cell growth and differentiation. Differentiation of endoderm-derived endocrine pancreas is mediated through PAX4 and PAX6. PAX4 may act as a PAX6 repressor due to the competition for binding sites and lower transactivation potential of Pax-4. PAX4 has the potential to function as a tumor suppressor in human melanoma.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



L02 cells were subjected to SDS PAGE followed by western blot with 21387-1-AP (PAX4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.