

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-TNFR1



Numéro de catalogue: 21574-1-AP

Phare

41 Publications

## Informations de base

<b>Numéro de catalogue:</b> 21574-1-AP	<b>Numéro d'acquisition GenBank:</b> BC010140	<b>Méthode de purification:</b> Purification par affinité contre l'antigène
<b>Taille:</b> 150ul , Concentration: 300 µg/ml by Nanodrop;	<b>Identification du gène (NCBI):</b> 7132	<b>Dilutions recommandées:</b> WB 1:500-1:1000 IHC 1:50-1:500 IF 1:50-1:500
<b>Hôte:</b> Lapin	<b>Nom complet:</b> tumor necrosis factor receptor superfamily, member 1A	
<b>Isotype:</b> IgG	<b>MW calculé:</b> 455 aa, 50 kDa	
<b>Immunogen Catalog Number:</b> AG16112	<b>MW observés:</b> 50 kDa	

## Applications

**Applications testées:**  
FC, IF, IHC, WB, ELISA

**Demandes citées:**  
CoIP, IF, IHC, IP, WB

**Spécificité de l'espèce:**  
Humain, souris

**Espèces citées:**  
Humain, rat, souris

**Contrôles positifs:**

**WB :** cellules Raji, cellules HeLa, cellules HL-60, tissu cérébral humain

**IHC :** tissu cérébral humain, tissu de cancer du sein humain

**IF :** tissu cérébral de souris,

**Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9,0; (\*) À défaut, le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6,0.**

## Informations générales

Tumor necrosis factor (TNF) is a multifunctional cytokine that plays a key role in regulating inflammation, immune functions, host defense, and apoptosis (PMID: 16407280). TNF exists in soluble and membrane-bound forms. TNF signals through two distinct cell surface receptors, TNFR1 (TNFRSF1A, CD120a) and TNFR2 (TNFRSF1B, CD120b). Whereas TNFR1 is widely expressed, expression of TNFR2 is limited to cells of the immune system, endothelial cells, and nerve cells (PMID: 22053109). TNFR1, which contains a death domain (DD) within its intracytoplasmic region, is thought to be the key receptor for TNF signaling (PMID: 16407280). This receptor can activate NF-kappaB, mediate apoptosis, and function as a regulator of inflammation. Antiapoptotic protein BCL2-associated athanogene 4 (BAG4/SODD) and adaptor proteins TRADD and TRAF2 have been shown to interact with this receptor, and thus play regulatory roles in the signal transduction mediated by the receptor.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Sisi Lei	36172180	Front Pharmacol	WB
Manli Wang	36099882	Cancer Cell	IF
Xian Wang	36147345	Front Pharmacol	WB,IHC,IF

## Stockage

**Stockage:**

Stocker à -20°C. Stable pendant un an après l'expédition.

**Tampon de stockage:**

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Données de validation sélectionnées




Raji cells were subjected to SDS PAGE followed by western blot with 21574-1-AP (TNFR1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.




$1 \times 10^6$  Raji cells were stained with 0.2ug TNFR1 antibody (21574-1-AP, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.




Immunohistochemical analysis of paraffin-embedded human brain tissue slide using 21574-1-AP (TNFR1 antibody) at dilution of 1:200 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human brain tissue slide using 21574-1-AP (TNFR1 antibody) at dilution of 1:200 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using TNFR1 antibody (21574-1-AP) at dilution of 1:200 and Coralite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using TNFR1 antibody (21574-1-AP) at dilution of 1:200 and Coralite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).