

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZNF132



Numéro de catalogue: 21681-1-AP

Informations de base

Numéro de catalogue:

21681-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop and 267 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG16142

Numéro d'acquisition GenBank:

BC109107

Identification du gène (NCBI):

7691

Nom complet:

zinc finger protein 132

MW calculé

706 aa, 81 kDa

MW observés:

81 and 68 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, souris

Contrôles positifs:

WB : tissu cardiaque de souris, tissu testiculaire humain

Informations générales

ZNF132, also named as Zinc finger protein 132, is a 706 amino-acid protein, which contains 18 C2H2-type zinc fingers and 1 KRAB domain. ZNF132 belongs to the krueppel C2H2-type zinc-finger protein family and localizes in the nucleus. ZNF132 may be involved in transcriptional regulation. ZNF132 exists as two alternatively spliced isoforms, and is encoded by a gene that maps to human chromosome 19.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse heart tissue were subjected to SDS PAGE followed by western blot with 21681-1-AP (ZNF132 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



human testis tissue were subjected to SDS PAGE followed by western blot with 21681-1-AP (ZNF132 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.