

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-SOX2



Numéro de catalogue: 21965-1-AP

## Informations de base

Numéro de catalogue: 21965-1-AP	Numéro d'acquisition GenBank: BC013923	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul	Identification du gène (NCBI): 6657	
Hôte: Lapin	Nom complet: SRY (sex determining region Y)-box 2	
Isotype: IgG	MW calculé 317 aa, 34 kDa	
Immunogen Catalog Number: AG16794		

## Applications

Applications testées:  
ELISA

Spécificité de l'espèce:  
Humain

## Informations générales

Sox2, also known as SRY (sex determining region Y)-box 2, is a transcription factor essential for maintaining self-renewal of undifferentiated ES cells and is one of the key transcription factors used to reprogram mouse and human fibroblasts to a pluripotent state. Sox2 expressed in undifferentiated pluripotent stem cells and germ cells during development. Affinity purified rabbit anti-Sox2 antibody can be used to demonstrate pluripotency of ES and iPS cells. This antibody is a rabbit polyclonal antibody raised against residues near the C terminus of human SOX2.

## Stockage

Stockage:  
Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

