

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-HTR1B



Numéro de catalogue: 22189-1-AP

1 Publications

Informations de base

Numéro de catalogue: 22189-1-AP	Numéro d'acquisition GenBank: BC096207	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 3351	Dilutions recommandées: WB 1:500-1:1000
Hôte: Lapin	Nom complet: 5-hydroxytryptamine (serotonin) receptor 1B	
Isotype: IgG	MW calculé 390 aa, 44 kDa	
Immunogen Catalog Number: AG16753	MW observés: 40-44 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : tissu placentaire humain, tissu placentaire de souris
Demandes citées: WB	
Spécificité de l'espèce: Humain, souris	
Espèces citées: souris	

Informations générales

HTR1B, also known as 5-HT1B, 5-HT1DB, HTR1D2, HTR1DB, and S12, belongs to the G-protein coupled receptor 1 family. HTR1B is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. The activity of HTR1B is mediated by G proteins that inhibit adenylate cyclase activity.

Publications notables

Autrice	Pubmed ID	Journal	Application
Huiduo Guo	34695613	Ecotoxicol Environ Saf	WB

Stockage

Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



human placenta tissue were subjected to SDS PAGE followed by western blot with 22189-1-AP (HTR1B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.