

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-AMOTL2



Numéro de catalogue: 23351-1-AP

Phare

5 Publications

Informations de base

Numéro de catalogue:

23351-1-AP

Taille:

150ul, Concentration: 700 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG19953

Numéro d'acquisition GenBank:

BC011454

Identification du gène (NCBI):

51421

Nom complet:

angiotensin like 2

MW calculé

779 aa, 86 kDa

MW observés:

100 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

Applications

Applications testées:

WB, ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain

Espèces citées:

Humain

Contrôles positifs:

WB : cellules HEK-293,

Informations générales

Angiotensin-like 2 (AMOTL2) is a member of the motin family of angiotensin-binding proteins. The motin family, also known as AMOTs, consists of three members: AMOT, AMOT-like 1 (AMOTL1) and AMOTL2. Members of the motin family are a type of adaptor proteins mainly distributed in the cytomembrane, cytoplasm or nucleus, and have a higher expression in the endothelial cells of capillaries and in larger vessels of the placenta. (PMID: 34036399, PMID: 25080976)

Publications notables

Autrice	Pubmed ID	Journal	Application
Diala Kantar	33707576	Sci Rep	WB
Patrizia Brunner	31800378	Mol Biol Cell	WB
Zengya Guo	32942233	Mol Ther Nucleic Acids	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HEK-293 cells were subjected to SDS PAGE followed by western blot with 23351-1-AP (AMOTL2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.